

# Spectra-H lux powered by EASTCON®

## Medium performance:

Radiation total	60-90 mW
Radiation Lase	0,4-6,25 mW
Radiation Infrared	30-90 mW
Red LEDs	2-10 mW

## Wavelengths:

Laser	905 nm
Infrared	875 nm
Redlight	635 nm

## General specifications:

Device type	Unit class II, application part type BF
Induction permanent magnet	25-45 mT
Output Laser	25 W
Laser category	I (IEC 60825-1)
Power supply	5 VDC, Max. 1,8A
AC/DC charger	100-240 V, 50-60 Hz
Operating mode	Permanent operation